

## Spencer 407/B

Transport chair with 4 wheels



Spencer 407/B is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

### Specific features

- Easy to clean PVC sheets
- Armrest and rear handles made of a single component, allowing weight reduction
- Robust and lightweight structure
- Footrest and armrests to improve the stability and patient comfort
- Large rear wheels to reduce vibration during patient transfer
- Swiveling front wheels with brake for rotation and axis
- Quick and intuitive closure system
- Safety system against accidental closure

### Reference standard

NA

### Standard equipment

ST10408A | 3 belts for chair 406/407

### Optional accessories

ST42705E | 10g wall fastener for chair 407

### Technical data

Length with handles closed	810 ± 10 mm
Length with handles opened	880 ± 10 mm
Width	475 ± 10 mm
Height opened	970 ± 10 mm
Height folded	800 ± 10 mm
Thickness folded	250 ± 10 mm
Rear wheels	∅ 200 mm
Front wheels	∅ 100 mm
Materials	Steel/Al/PVC
Weight	12,5 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 407B CHROME  
W/ARMRESTS EVACUATION CHAIR  
4 WHEELS



**ST40407A**

**CND Classification** V9099

**Registration number** 93563

Rev.0 (05/07/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)