

HIRVIN

Carrying sheet



Hirvin, is a device intended to transfer patients from the stretcher to the bed or for transport from the rescue place to the stretcher for patients who do not have spinal injuries.

Specific features

- Built with fire-resistant materials
- Easy sanitation thanks to the non-absorbent materials
- Nylon stitching to ensure high resistance
- The width of the tape of the handles, ensure a secure and safe grip
- Require minimal storage space

Technical data

Length	1850 ± 20 mm
Width	760 ± 5 mm
Dimensions wrapped ⁽¹⁾	150x80x760 mm
N° handles	8
Tapes width	50 mm
Handles lenght	170 mm
Loading capacity	150 kg
Materials	Nylon, PVC

⁽¹⁾ The measure refers to the wrapped device. It is just an example, because the device can be folded or wrapped in many ways.

Accessories

- | | |
|----------|-------------------------------|
| CB03000A | Carrying bag orange w/handles |
|----------|-------------------------------|

Class I Med.Device according to Dir. 93/42/CEE

HIRVIN TRANSFER SHEET 8
HANDLES
QM40200A



CND Classification	V0899
Registration number	170944

Reference standards

UNI EN 1865-1

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it