

HP HIRVIN PLUS

High resistance carrying sheet



Hirvin plus, is a device intended to transfer patients from the stretcher to the bed or for transport from the rescue place to the stretcher for patients who do not have spinal injuries.

Specific features

- Built with fire-resistant materials
- Easy sanitation thanks to the non-absorbent materials
- Nylon stitching to ensure high resistance
- The width of the tape of the handles, ensure a secure and safe grip
- Require minimal storage space
- Dedicated compartment to store a pillow or the supplied carrying bag

Accessories supplied

- CB03000A Carrying bag orange w/handles

Reference standard

UNI EN 1865-1

Technical data

Length	1850 ± 20 mm
Width	760 ± 5 mm
Dimensions wrapped ⁽¹⁾	150x80x760 mm
N° handles	8
Tapes width	50 mm
Handles lenght	170 mm
Loading capacity	150 kg
Materials	Nylon, PVC
Weight	1,5 ± 0,2 kg

⁽¹⁾ The measure refers to the wrapped device. It is just an example, because the device can be folded or wrapped in many ways.

Class I Med.Device according to Dir. 93/42/CEE

HP HIRVIN PLUS YELLOW
W/BAG
QM40501A



CND Classification	V0899
Registration number	170954

Rev.3 (09/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it