

FP

Rear fastener for Carrera Pro



System to be used to fix the stretchers on the floor of the medical vehicle. The device is to be used coupled with the corresponding front fixing. This version is specific for stretchers Carrera Pro.

Specific Characteristics

Dimensions reduced to the minimum
Coupling mechanism simple and intuitive

Technical data

Dimensions ^(hook up)	115x130x80 ± 2mm
Dimensioni ^(hook down)	115x170x37 ± 2mm
Materials	Al, Steel, cast iron
Weight	1,16 ± 0,1 kg

Class I Med.Device according to Dir. 93/42/CEE

F.P. REAR FASTENER FOR CARRERA
PRO

ST42107A

CND Classification



V9099

Registration number

789127

Rev.0 (25/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it