

QMX777

Anatomic Mattress



Mattress intended for use with stretchers, in order to increase patient comfort by absorbing the stresses to which it may be subjected during the transfer stages.

Specific features

- Smooth and easy to clean surface
- Tape with strap elements, used to fix the mattress to the stretcher board
- The shape, the internal divisions and the filling material, allow excellent adaptation to the patient's body
- Compatible with boards type S and T

Accessories

DR60851K DRX 850 BOX 100 PCS - TNT SHEET

Technical data

Dimensionis	1770x430 ± 10 mm
Central part thickness	50 mm
Maximum thickness	90 mm
Materials	PU exp, PVC,
Weight	2,3 ± 0,1 kg

Class I MD. According to Dir. 93/42/CEE

QMX 777 BLACK ANATOMIC
MATTRESS IN PVC
ST7000A



CND Classification	V9099
Registration number	93126

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it