

QMX 750

Mattress with handles



A mattress that can be used also for transfer or patient thanks to the integrated handles.

Specific features

Smooth surface easy to clean

Zip on the upper side to allow the extraction of the padding making cleaning procedures easier

Wide handles for a comfortable grip

8 handles allow the use by more operators reducing their effort

Technical data

Dimensions	1850x520 ± 10 mm
Thickness	50 mm
Materials	PU exp, PVC
Weight	2,55 ± 0,1 kg

Class I MD. According to Dir. 93/42/CEE

QMX 750 MATTRESS WITH
HANDLES

QM50150A

CND Classification V9099

Registration number 93009



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it