

QMX 750

Mattress with handles



A mattress that can be used also for transfer or patient thanks to the integrated handles.

Specific features

- Smooth surface easy to clean
- Zip on the upper side to allow the extraction of the padding making cleaning procedures easier
- Wide handles for a comfortable grip
- 8 handles allow the use by more operators reducing their effort

Technical data

Dimensions	1850x520 ± 10 mm
Thickness	50 mm
Materials	PU exp, PVC
Loading capacity <small>(distributed on all 8 handles)</small>	120 kg
Weight	2,55 ± 0,1 kg

Class I MD. According to Dir. 93/42/CEE

QMX 750 MATTRESS WITH
HANDLES

QM50150A

CND Classification V9099

Registration number 93009



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it