


QMX 02

Multilevel Pillow



Pillow with adjustable thickness, to be installed on the backrest of Spencer self-loading stretchers.

Specific Characteristics

Easy to clean
 The knob for thickness adjustment has an ideogram for an intuitive use
 Seamless integration with the backrest of the stretcher
 Low weight

Technical data

Dimensions	330x300 mm
Materials	PVC, PE
Weight	970 g
Thickness adjustable	from 35 to 75 mm

Class I Med.Device according to Dir. 93/42/CEE

MULTI LEVEL PVC PILLOW FOR
 DNA BOARDS
QM50002A



CND Classification	V0899
Registration number	171905

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it