

ZEN

Ergonomic cushion



Zen is an ergonomic cushion designed to provide support, to vary the head inclination from the patient floor, increase patient comfort and to contribute to hyperextend the head improving airway's opening.

Specific features

- Smooth surface easy to clean
- The special shape allows dual mode of application in order to offer different types of support
- The strap bands on the built-in tape, confers stability when placed on surfaces equipped with male strap
- The tape allow to fix the pillow to handle or to a tubular part of the patient board
- The raised edges allow a perfect coupling with Spencer mattresses
- Applied with the higher part in the direction of the base of the neck, allows hyperextension of the head increasing airway's opening

Accessories

DR60550K	WHITE DISPOSABLE COVER FOR ZEN 37.2X34CM (50PCS)
----------	---

Technical data

Dimensions	235x300x130 ± 10 mm
Materials	PU, Nylon, polyester
Weight	680 ± 20 g

Class I Med.Device according to Dir. 93/42/CEE

ZEN ANATOMIC PILLOW YELLOW
QM50000A



CND Classification	V0899
Registration number	170966

Rev.2 (18/07/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it