

## SXV 20G

### Vertical fastener for Pick Up stretcher



Fastener studied for a safe fixing of Spencer SX Pick Up stretchers.

#### Specific Characteristics

Extremely robust and lightweight  
 Easy to use and install  
 Designed to withstand accelerations up to 20 g

#### Technical data

Dimensions <sup>(long piece)</sup>	440x25x95 ± 5 mm
Dimensions <sup>(short piece)</sup>	180x60x60 ± 5 mm
Material	Al, Nylon, brass
Weight <sup>(overall)</sup>	790 ± 10 g

#### Compliant with

EN 1789

#### Class I Med.Device according to Dir. 93/42/CEE

SXV VERTICAL WALL BRACKET 20  
 FOR PICK UP STRETCHER

**ST05018C**

**CND Classification** V9099

**Registration number** 91914



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)