

## STX 40

### Shell carry bag



Bag for transport and storage of basket stretchers Spencer Shell.

#### Specific Characteristics

PVC fabric easy to clean

Handles coated with soft plastic material for easy transport

Longitudinal belts with ring terminal allow wall hanging

Upper strap closure

#### Technical data

Dimensions	2150x610x200 ± 10 mm
N° of handles	6
Material	PVC
Weight	2,35 ± 0,05 kg

STX 40 BAG FOR SHELL

**ST04040A**

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)