

F 400

20 G fastener



The fastener is intended for installation inside a medical vehicle. The device is compatible with Spencer evacuation chairs of Profile and Weld series.

Specific Characteristics

- Compact and robust
- Foot operated simple and intuitive mechanism
- Fixing belt with integrated reel
- Quick release metal buckle

Technical data

Base dimensions	401x72x80 mm
Upper part dimensions	480x30x60 mm
Materials	Steel, Nylon
Weight	1,9 kg

Compliant with

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

F 400 FIXATION SYSTEM 20 G FOR CHAIRS

ST36451B

Classification CND V9099

Registration number 93265



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it