

DNA STRAP

Belt with integrated re-winding system



Adjustable belt to stabilize the patient during transport. The reel body is equipped with an interlocking structure which perfectly fits into the patient boards of Spencer stretchers.

Specific features

- Semi automatic reel to increase rapidity during intervention
- Quick coupling system with metal buckle

Technical data ¹

Length	Max 1360 mm
Belt width	50 mm
Materials	Nylon, PP, Steel
Weight	380 ± 10 g

¹ Indicated length is subject to tolerance of ± 30 mm

Class I Med.Device according to Dir. 93/42/CEE

DNA STRAP BELT W/INTEGRATED
RE-WINDING SYSTEM

ST00497B

CND Classification



V9099

Registration number

96175

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it