

DNA STRAP THORAX BELT

Belt with integrated re-winding system



Adjustable thorax belt to stabilize the patient during transport. The reel body is equipped with an interlocking structure which perfectly fits into the patient boards of Spencer stretchers.

Specific features

- Semi automatic reel to increase rapidity during intervention
- Quick coupling system with metal buckle

Technical data ¹

Central belt length	Max 1360 mm
Shoulder belt length	Max 1150 mm
Belt width	50 mm
Materials	Nylon, PP, Steel
Weight	920 ± 10 g

¹ Indicated length is subject to tolerance of ± 30 mm

Class I Med.Device according to Dir. 93/42/CEE

DNA STRAP THORACIC BELT W/RE-WINDING SYSTEM
ST00498B



CND Classification	V9099
Registration number	96176

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it