

# STX 597

## One piece belt



Adjustable belt to stabilize the patient during transport procedures.


### Specific features

Nylon hook to ensure maximum safety  
Suitable for spinal boards, pick up stretchers or other type of fastening

### Technical data <sup>1</sup>

	Min (mm)	Max (mm)
Length L	710	1600
Belt width	50 mm	
Materias	Nylon, PP	
Belt weight	130 ± 10 g	

<sup>1</sup> The indicated lengths are subject to tolerances of ± 30 mm



Class I Med.Device according to Dir. 93/42/CEE

STX 597 – ONE PIECE BETLS

ST00597A

CND Classification



V9099

Registration number

95879

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)