

STX 591

One piece belt



Belt with metal buckle used to ensure a proper containment of the patient. The device is suitable for use with vacuum mattresses and pick-up stretchers.


Specific features

- Extremely versatile and easy to use thanks to the single metal buckle with quick release system
- Characterized by a wide adjustment range
- Metal buckle to ensure maximum safety

Technical data ¹

	Min (mm)	Max (mm)
Length L	350	1850
Belt width	50 mm	
Materials	Steel, Nylon	
Weight	250 g	

¹ The measures are the maximum reachable. The effective length can vary depending on the device and the knot done to fix the belt to the stretcher



Class I Med.Device according to Dir. 93/42/CEE

MOD 591 RESTRAINT 1 PIECE

METAL

ST00591A

CND Classification

Registration number



V9099

96040

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it