

STX 702

Metal buckle belt



Couple of adjustable belts intended to stabilize patient during transport. The belt must be anchored to the frame of the stretcher by means of a loop knot.


Specific features

- Metal buckle to provide more safety
- Reflective insert inside the fabric to increase visibility
- Suitable for self loading and basket stretchers

Technical data ¹

	Min (mm)	Max (mm)
Length L	680	1500
Tape width	50 mm	
Materials	Steel, Nylon	
Weight of the pair	460 g	

¹ The indicated sizes are in reference to the maximum dimensions. The useful length can be different, due to the device ore to the knot used to attach the belt to the frame of the stretcher.



Class I Med.Device according to Dir. 93/42/CEE

STX 702 BLACK BELT TWO PIECE,
METAL BUCKLE, REFLEX
ST70002A



CND Classification	V9099
Registration number	96045

Rev.3 (17/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it