

STX 598


Two piece belt



Adjustable belt to stabilize the patient during transport procedures. It must be anchored to the frame of the stretcher. The fastening must be done by means of a loop knot .

Specific features

Nylon hook to ensure maximum safety
Suitable for self loading stretchers, basket stretchers and spinal boards



Technical data ¹

	Min (mm)	Max (mm)
Length L	710	1500
Available length after knot has been done ²	460	1250
Belt width	50 mm	
Materias	Nylon, PP	
Belt weight	130 ± 10 g	

¹ The indicated lengths are subject to tolerances of ± 30 mm

² Approximate value. It varies depending on which device the belt is applied

Class I Med.Device according to Dir. 93/42/CEE

MOD 598 RESTRAINT BELT 2 PCS.

DERLIN

ST00598A

CND Classification

Registration number



V9099

96090

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it