

STX 592

Two pieces belt



Adjustable belt to stabilize the patient during transport procedures. It must be anchored to the frame of the stretcher. The fastening must be done by means of a loop knot.

Specific features


- Metal hook to ensure maximum safety
- Suitable for self loading stretchers, basket stretchers and spinal boards

Technical data ¹

| | Min (mm) | Max (mm) |
|--|------------|----------|
| Length L | 660 | 1580 |
| Available length after knot has been done ² | 410 | 1330 |
| Belt width | 50 mm | |
| Materials | Nylon, PP | |
| Belt weight | 250 ± 10 g | |

¹ The indicated lengths are subject to tolerances of ± 30 mm

² Approximate value. It varies depending on which device the belt is applied



Class I Med.Device according to Dir. 93/42/CEE

STX 592 - TWO PIECE METAL BELT

ST00592A

CND Classification



V9099

Registration number

96104

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it