


STX580 - SPEED BASKET

Basket carrier



Speed Basket is a supporting system designed to facilitate the movement of Spencer basket stretchers on the ground.

Caratteristiche specifiche

Robust and lightweight structure

Large diameter wheel to make transport easier

The positioning of the basket stretcher is quick and easy thanks to the 4 rapid locking systems

Technical data

Height	830 ± 10 mm
Width	670 ± 10 mm
Length	840 ± 10 mm
Materials	Al, steel, Nylon
Loading capacity	160 kg

Class I med.dev. according to dir. 93/42/CEE

STX580 - BASKET STRETCHER

CARRIER

ST04600B

CND Classification



V9099

Rev.0 (22/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it