

K BELT 3

Belts one piece with bag




Adjustable belt to stabilize the patient during transport procedures. This kit comes with a carry bag.

Specific features

- Nylon hook to ensure maximum safety
- Suitable for spinal boards, pick up stretchers or other type of fastening
- Transport bag for a neater storage

Technical data ¹

	Min (mm)	Max (mm)
Length L	710	1600
Belt width	50 mm	
Materials	Nylon, PP	
Single belt weight	130 ± 10 g	
Kit weight	450 ± 20g	



¹ The indicated lengths are subject to tolerances of ± 30 mm

Class I Med.Device according to Dir. 93/42/CEE

K BELT 3 - THREE ONE PIECE BELTS
STX 597

ST00593A

CND Classification



Registration number

V9099

98851

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it