

K BELT 1

Belts two piece with bag



Adjustable belt to stabilize the patient during transport procedures. It must be anchored to the frame of the stretcher. The fastening must be done by means of a loop knot. This kit comes with a carry bag.

Specific features


- Nylon hook to ensure maximum safety
- Suitable for self loading stretchers, basket stretchers and spinal boards
- Transport bag for a neater storage

Technical data ¹

	Min (mm)	Max (mm)
Length L	710	1500
Available length after the knot has been done ²	460	1250
Belt width	50 mm	
Materials	Nylon, PP	
Single belt weight	130 ± 10 g	
Kit weight	445 ± 20g	

¹ The indicated lengths are subject to tolerances of ± 30 mm

² Approximate value. It varies depending on which device the belt is applied



Class I Med.Device according to Dir. 93/42/CEE

K BELT 1 - THREE TWO PIECE BELTS

STX 598

ST00594A

CND Classification

Registration number



V9099

98849

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it