

JACOB

Anti-shock trousers



Device intended for use on patients with pelvic injuries.

Specific features

- Equipped with straps for easy application and adjustment
- 3 separate air chambers for independent adjustment of the inflation pressure
- The pump allows to inflate all chambers at the same time
- Large inguinal opening to allow access to the femoral vessels
- The transport bag can easily store the device together with its accessories
- 4 lifting handles for easy handling of the patient
- Calibrated overpressure valves to prevent an inflating of chambers to dangerous levels

Technical data

Dimensions	A: 530 ± 20mm B: 340 ± 20mm C: 780 ± 20mm
Bag dimensions	550x330x160 ± 10mm
Trousers weight	2,35 ± 0,3 kg
Pump weight	0,68 ± 0,3 kg
Manometer scale	0-300 mmHg
Pump tube length ^{Main tube}	120 ± 2cm
Materials	Nylon, PU, PVC



Standard di riferimento

NA

Class I Med.Device according to Dir. 93/42/CEE
JACOB ANTISHOCKS PANTS ADULT



CND Classification V9099
Registration number 195452

Rev.1 (10/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it