

CROSS CHAIR

Trolley 10G for Stretcher/Chair



Trolley designed to be used with the Stretcher/Chair Cross Chair, for easier handling, loading and fixing inside ambulances. 10G certified if used with Spencer fastener R-Max or S-MAX.

Specific Characteristics

- Robust stainless steel frame
- Wheels $\varnothing 200$ to increase sliding action
- Parking brakes on rear wheels
- Possibility to place in intermediate height to facilitate loading of Cross Chair
- Cross Chair fixing system easy to operate

Accessories

CR00080C	Cross Chair Stretcher/Chair
ST42707B	Fastener S-Max <small>(Mandatory for 10g versions)</small>

Technical data

Lenght	1980 \pm 10 mm
Width	570 \pm 5 mm
Sliding profile height (Raised trolley)	850/900 mm
Sliding profile height (Trolley in intermediate height)	600/680 mm
Sliding profile height (Lowered trolley)	280 mm
Wheels	\varnothing 200 mm
Weight	35,5 \pm 0,5 kg
Loading capacity ¹	170 kg

¹ Considering patient weight, equipment and accessories

Class I Med.Device according to Dir. 93/42/CEE

CROSS CHAIR TROLLEY 10G

CR00082C

CND Classification



V0899

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it