

## CROSS CHAIR TROLLEY

Trolley for Stretcher/Chair



Trolley designed to be used with the Stretcher/Chair Cross Chair, for easier handling, loading and fixing inside ambulances. It must be used with the dedicated Spencer fixing systems .

### Specific Characteristics

- Robust stainless steel frame
- Wheels  $\varnothing 200$  to increase sliding action
- Parking brakes on rear wheels
- Possibility to place in intermediate height to facilitate loading of Cross Chair
- Cross Chair fixing system easy to operate

### Technical data

Lenght	1980 $\pm$ 10 mm
Width	570 $\pm$ 5 mm
Sliding profile height (Raised trolley)	850/900 mm
Sliding profile height (Trolley in intermediate height)	600/680 mm
Sliding profile height (Lowered trolley)	280 mm
Wheels	$\varnothing$ 200 mm
Weight	35,5 $\pm$ 0,5 kg
Loading capacity <sup>1</sup>	170 kg

<sup>1</sup> Considering patient weight, equipment and accessories

### Accessories

CR00080C	Cross Chair Stretcher/Chair
ST42100A	Posterior fastner (FP)
ST42200A	Anterior fastner (FA)

Class I Med.Device according to Dir. 93/42/CEE

CROSS CHAIR TROLLEY

CR00081C

CND Classification



V0899

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)