

PIN STRAPS

Belt system with hooks



Retention system suitable for use with all spinal boards equipped with at least 5 pins and 5 handles for each side.

Specific features

Fastening system with carabiner for quick and easy fastening to spinal boards equipped with pins
 Strap closure system to allow the use also with spinal boards without pins
 4 transversal belts 3 of which can slide for a better adaptation to device in use

Technical data ¹

Lengths	Min (mm)	Max (mm)
Length of side belts	580 ± 30	970 ± 30
System vertical length	1160 ± 50 mm	
Materials	Nylon/PP/Al	
Weight	1500 ± 20 g	

¹ Lengths are referred to the installed device. The real length can vary depending on the device in use with the product..

Class I Med.Device according to Dir. 93/42/CEE

PIN STRAPS- SYSTEM OF BELTS
 WITH HOOKS

ST02015A

CND Classification V9099

Registration number 99840



Rev.1 (04/12/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it