


STX 90

Telescopic headrest



Anatomic pillow with telescopic support, to be used with DNA patient floors of Spencer stretchers. The device provides neck support, increases patient comfort and improves head alignment during patient transport allowing also to hyperextend the head.

Specific Characteristics

- Smooth surface easy to clean
- Its specific shape allows to apply the pillow in two ways, providing different types of support
- Continuous adjustment for maximum adaptability
- The raised edges allow a perfect fit with Spencer mattresses
- If installed with the high part of the pillow oriented toward the base of the neck, allows to hyperextend the head improving airway's opening

Technical data

Dimensions	480x300x230 ± 10 mm
Materials	PU, Steel, Nylon,
Weight	1050 ± 20 g

Accessories

DR60550K	ZEN PILLOWCASES BOX 50 PCS - ERGONOMIC CUSHION
----------	---

Class I Med.Device according to Dir. 93/42/CEE

STX 90 TELESCOPIC HEAD REST
FOR STRETCHERS

ST00491A

CND Classification



V0899

Rev.1 (18/07/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it