

REK-OX

Manual switcher



Manual exchanger intended for use with Spencer oxygen delivery systems. The device is designed for use with Rek-Ox systems and allows manual selection of the supply oxygen tank simply by turning the knob.

Specific features

- Maximum integration with the distribution system
- The selection is intuitive thanks to the big sized identification numbers and to the colour coded positions
- Compact design achieved thanks to a very small pneumatic valve
- High structural strength thanks to the one-piece aluminum body

Technical data

Dimensions	90x70x65 mm
Materials	Al, PP, ABS
Input and output connections	M 1/8
Working pressure	400 ₀ ⁺¹⁰⁰ kPa
Weight	0,5 kg

Class IIa Med. Dev. According to Dir. 93/42/CEE

REK-OX SWITCH
BOTTLE.1/BOTTLE.2
OX07202



CND Classification	V9099
Registration number	1271659

Rev.0 (21/12/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it