

HUX

Fastener for stretcher



Fastener for stretchers, created to provide an anchorage inside medical vehicles.

Specific features

The linearity of the structure avoids unnecessary bulk when the stretcher is not loaded on the ambulance

The release of the stretcher happens easily, simply operating on the small release mechanism on the rear fastener

The device acts simultaneously both on frontal and rear part of the stretcher and, thanks to the integral structure, allows a simple installation without alignment difficulties between the frontal and rear part of the fastener

The weight is minimized thanks to researches of advanced materials and in-depth studies about constructive geometries

The coupling between the rear fastener and the main body, allows easy installation thanks to the possibility of adjustments

Technical data

Dimensions	1625x190x80 ± 20mm
Metal profile thickness	30 mm
Materials	Al , Nylon, PE, Steel
Weight	5,6 ± 0,2 kg

Class I Med.Device according to Dir. 93/42/CEE

HUX INTEGRAL FASTENER FOR STRETCHERS

ST42712

CND Classification



V9099

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it