



Re-organise your work while at the same time protecting and transporting your most sophisticated products and instruments. K Series is a new line of bags using the Matrix body created by Spencer for the emergencies of today and of the future. K Series: a collection of bags and novelties which have been studied to be practical, functional and resistant to all atmospheric conditions as well as the wear and tear of time. For a rapid identification, each bag has a different colour code.

External covering: scratch proof and washable expanded polyethylene

Internal: nylon taffetà, elasticized loops, transparent PVC

Closure: zip complete with zip pull

Dedicated to the Mark backpack and Smart bag. Suitable for bags, backpacks and cases.

TECHNICAL DATA

Dimensions (mm):	295x65xh180
Capacity (L):	3
Weight (g):	300
External material:	PU
Heat formed:	Yes
Internal material Taffetà	
Nylon/Cordura:	Yes
Zipper closure:	Yes
Zip pull:	Yes
Compartments / internal pockets:	2
Elastic loops:	60
Detachable internal flap:	Yes
Pharmaceuticals and subministration:	Yes

K 6 red CB04143B