

E-MAX

10G fastener for EVA



Fastener intended for installation inside an ambulance, to be used with EVA chairs. The device should be fixed on wall, according to the vehicle structure and equipment.

Specific Characteristics

The polyethylene base, dampens vibrations and makes cleaning procedures easier

The locking of the chair take place simply by placing the chair on the base and pushing it against the fastener

Characterized by high resistant and slim body

Technical data

Dimensions	425x320x155 ± 5mm
Materials	PE, Nylon, Al
Weight	3,5 ± 0,2 kg

Reference standards

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

E-MAX WALL BRACKET 10G TESTED FOR EVA CHAIR

SK30001E



CND Classification V9099

Registration number 445722

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it