

Res-Q-Matt plus

Vacuum mattress



Vacuum mattresses are devices for the immobilization of trauma patients. The vacuum system allows different levels of immobilization. The more rigid helps the right positioning of the patient; softer ones allow to mitigate the low-frequency vibrations typically found on emergency vehicles during transport.

Specific characteristics

Extremely resistant fabric
Integrated handles allow a firm grip

Reference standard

NA

Optionals

QM22120A	Aluminum Pump
QM22125A	Aluminum double effect pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit

Technical data

Length	2140 ± 20mm
Width	940 ± 10mm
Filling material	Polystyrene
N° of handles	8
Materials	Al, PVC, Nylon
Weight	5 ± 0,5 kg
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-MATT PLUS POLYSTYRENE-
PROFESS. VAC. MATTR.

QM22101A

CND Classification



V0804

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it