

## GMX 223

### Surgical straight scissors



Surgical straight scissors pointed end.

#### Specific features

Sterilizable

The build quality allows a precise cut

#### Technical data

Overall length	140 mm
Blade length	45 mm
Weight	38 g
Material	Stainless steel

Class I Med.Device according to Dir. 93/42/CEE

SURGICAL STRAIGHT SCISS

POINTED END 14.5 cm

GM00223A



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)