

NXT MAX P

10g fastener/charger

For NXT lung ventilators



10g fastener and battery charger with quick release system for use with Spencer NXT electronic lung ventilators

Specific features

Its shape is designed to allow easy and quick placing of the lung ventilator on the fastener

The quick release button on the front side is easy to operate

Thanks to the perfect integration between the fastener and the lung ventilator, the overall dimensions remain almost the same as those of the lung ventilator alone, increasing its height by only a couple of centimeters

Easy to install inside the emergency vehicle, is fastened using the base

Compatible with 170NXT and 190NXT lung ventilators

The battery charging starts simply by placing the lung ventilator on the fastener

Reference standard

EN 1789:2007 + A2:2014

Technical data

Dimensions (DxWxH)	160x160x95 ± 5mm
Materials	HDPE, Nylon, Steel
Approx Weight	430 ± 10g
Electrical connection	12V Cigarette lighter

Class I Med.Device according to Dir. 93/42/CEE

NXT MAX P - 10G QUICK RELEASE
FASTENER/CHARGER

EV50121B



Rev.0 (12/02/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza(PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it