

NXT MAX WALL

10g wall fastener

For NXT lung ventilators



10g wall fastener for use with Spencer NXT lung ventilators

Specific features

Its shape is designed to allow easy and quick placing of the lung ventilator on the fastener

Thanks to the perfect integration between the fastener and the lung ventilator, the overall dimensions remain almost the same as those of the lung ventilator alone, increasing its height by only a couple of centimeters

Easy to install inside the emergency vehicle, is fastened using the back side, ideal solution for space optimization

Compatible with 118NXT, 170NXT, 190NXT lung ventilators

Reference standard

EN 1789:2007 + A2:2014

Technical data

Dimensions (DxWxH)	160x160x95 ± 5mm
Materials	HDPE, Steel
Approx Weight	300 ± 10g

Class I Med.Device according to Dir. 93/42/CEE

NXT MAX WALL - 10G WALL

FASTENER

EV50122B



Rev.0 (12/02/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza(PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it